

PPE

Breathable PPE Suit Hooded Protective Coverall with Shoe Cover for Corona Warriors



SALIENT FEATURES

- **Light Weight** - Complete suit including shoe cover, weighing only 300 gms approx.
- **Unique barrier coating** on outside and smooth, flexible and comfortable nano coating on inside, which allows **moisture vapour transmission & breathability**.
- Anti- Bacterial finish on the outside which disallows bacteria.
- Closely Woven, Compact Polyester Fibre of High Strength.
- **Seams are sealed** with sealing tape using a high frequency automatic machine
- *Design, **Developed and manufactured** in association with **DMSRDE, Kanpur (Wing of DRDO)***

- High Level of protection for Healthcare Professionals and frontline response personnel in most impacted regions of Corona Virus
- Blood Penetration test as per ISO 16603:2004 with highest rating
- Unique certificate code (UCC Code)
BPT000104XXXX
- Can be used after disinfecting for upto 5 times
- Comes in **various colours** for nurses, police personnel and safai karmi.

Disinfecting Procedure

First method

1. Disinfect by dipping in a cold bleach solution.
5 teaspoons of home bleach in 5 litres water
(5 / 6 % Chlorine)
Do not wriggle wash. Plain dipping for 2/3 minutes
2. Rinse with normal water
3. Air dry


Second method

1. Fogging of used coverall with 11% Hydrogen Peroxide solution in a sealed room.

Note: Use upto 5 times. Do not use IPA based disinfectants.



**HILTEX
PPE
COVERALL**



**Breathable & Water
Proof Coating**



Hood Back View



Seam Sealed Finishing



Long Shoe Cover



Taped Front Flap



Best Quality PPE Kit developed by DMSRDE (DRDO Kanpur)

5 times reusable + and most economical, only Rs.200/- per use


Comparison between Non-Woven Suit and Our Suit

Parameters	Non-Woven Laminated	Our Polyester woven fabric (reusable)	Remarks
Weight of PPE Kit	500 gms	300 gms	Our suit is light weight by 40%
Strength	100 N(gets torn easily)	300 N	strength is higher by 300%
Breathability	Zero	Dissipates body heat	Provides Comfort
Moisture Vapour Transmission	No	Yes	Provides Comfort
How many times can use	One time	Upto 5 times after disinfecting	Cost per use is as low as Rs. 200
Biodegradable	Has to be burned	Recyclable	Eco-friendly
Anti-Bacterial	NO	Yes	Provides higher safety

Specification for fabric for making Anti Viral Nano Coated Fabric for making Bio / Medical Suit and Shoe Cover

Sl No.	Testing Parameter	Requirement	Method of Test
1	Nature of Material	Polyester /Nylon	IS: 667
2	Mass of Fabric	60Gsm to 90 Gsm	IS: 7016 Pt. 1
3	Tensile Strength N/5	Warp > 450	IS: 7016 Pt. 2
		Weft > 200	
4	Tearing Strength	Warp > 15 N	IS: 7016 Pt. III
		Weft > 15 N	
5	Water Repellancy	>80	IS: 390
6	Water proofness (Coated Side) * Note : observation in respect of wetting etc. with in 5 mm of inside of gasket will be ignored.	No percolation of water through the fabric or wetting of the outer surface*	IS: 7016 Pt. VII
7	Water Vapour Transmission (Breathability) Option 'B', Temp-230C, Humidity in cup - 100%, Humidity in chamber- 50% and Air flow - 2.5 m/sec	2 mg/cm ² /hour (max.)	As per ASTM-96-95

FAX / SPEEDPOST / REGISTERED/BY HAND

Telephone	: 0512-2451759-78		No:- DRDO/DMSRDE/ACS/Coverall/01
Telegram	: LABDEV		Government of India, Ministry of Defence
Telex	: 325-385 DMSRD IN		Defence Research & Development Organization
Fax	: 0512-2453212		Defence Materials & Stores Research & Development Establishment
email : dmsrde@sanchar.net.in			DMSRDE P.O., GT Road, Kanpur - 208013
All correspondence to be addressed to the Director and this Estt. Reference be quoted invariably, while replying.			Date : 23.04.2020

To

M/s HILLTEX Industrial Fabrics Pvt. Ltd.
26/50, Birhana Road, Kanpur, India

Sub: **Issue of Unique Certification Code (BPT0001042020)**

Dear Sir

Your industry's developed '**Anti-Covid SUIT-I**' based on the polyester material under the technological support of DMSRDE, Kanpur is now tested for synthetic blood penetration test at Small Arm Factory (SAF), Kanpur vide our establishment Letter No. DRDO/DMSRDE/ACS/01 Dated 12/04/2020. The test result is indicating excellent level of protection against the penetration of Corona Virus and as the material is qualifying the requirement, a Unique Certification Code (**BPT0001042020**) is generated for the Suit (both fabric and sealed & stitched seam). This UCC is specifically allotted to this product (**Anti-Covid SUIT-I**) only.

This is for your kind information, please and your valuable support in this crucial phase of 'Corona Virus Pandemic' is highly acknowledged.

BR Das
23/04/2020
(Dr. Biswa Ranjan Das)

Sc 'D', DTT (DMSRDE), Kanpur

For Director
Scientist 'D' & Asst. Director
Govt. of India, MoD (DRDO)
DMSRDE, G.T. Road, Kanpur-208013

	 Small Arms Factory, Kalpi Road, Kanpur Ph. No. 0512-2295042-46, Fax No. 0512-2296229
	Report on Resistance of Materials used in Protective Clothing to Penetration using Synthetic Blood Test Report No SAF/QC(LG)/Lab/0001, Date: 17.04.2020
Sample details with Ref. no. and date	ANTI-COVID SUIT-I No. DRDO/DMSRDE/ACS/01, Dated 12/04/2020
Manufacture/Developer	DRDO, GT Road Kanpur-208013

- 1. Aim of the test:** To determine the resistance of protective clothing materials to penetration by blood and bodily fluids using Synthetic Blood
- 2. Equipment used for testing:** Blood Penetration Resistance Testing Equipment
- 3. Test Reference:** ISO 16603 - 2004 ; Unique Certification Code (UCC) - BPT0001042020
- 4. Materials used for testing:**

Synthetic Blood:

S. No.	Parameter	Value
1	Surface Tension	0.040 N/m
2	pH	7.22

5. Environmental Conditions

S. No.	Parameter	Value
1	Temperature	25°C
2	Relative Humidity	65%

6. Test Results:

S. No.	Sample Nomenclature	Sample Code	#Result as per ISO 16603-2004, Procedure C
1	Coverall fabric	A	5/6 (PASS)
2	Coverall Arm Sleeve	A	5/6 (PASS)
	Seam		
	B- Coverall Side Seam	A	N/A
	C- Coverall Leg Seam	A	6/6 (PASS)

#Samples have been tested as submitted and may not represent the batch/lot

Note: The procedures enumerated at Para-3 and 4 of the instructions of Ministry of Textiles, Government of India issued vide letter no. F. No. 8/4/2020-R&D dated 06.04.2020 must be strictly complied with regards to use of this certification.


Examiner


Authorized Signatory 17.4.2020